

Level 0
0 - <10 min

Level 1
10 - 30 min

Level 2
30 - 60 min


Level 3
60 - 120 min

Level 4
120 - 240 min

Level 5
240 - 480 min

Level 6
< 480 min

SHIELDskin XTREME* Sterile ORANGE NITRILE* 300 DI



- Category III PPE glove (PPE Directive 89/686/EEC)
- Complex Design - For mortal & irreversible risks
- Powder-free orange nitrile/polychloroprene glove
- Hand-specific
- 300 mm / 0.12 mm (EN 420:2003 + A1:2009)
- AQL 0.65 (EN 374-2:2014 Level 3)
- Viral penetration test (ISO 16604:2004 Procedure B & ASTM F1671-97b)
- Waterproof and for low chemical protection
- Tested for chemical permeation (EN 16523-1:2015 supersedes EN 374-3:2003)
- Typical particle levels: less than 2500 per cm² more or equal 0.5µm
- Type I hypersensitivity eliminated - Type IV hypersensitivity reduced

Mixed Solution
SPOR-KLENZ Solution

LEVEL 6
480 min

Mixed Solution
Perform sterile PAA ready to use

LEVEL 6
480 min

Mixed Solution
Perform sterile concentrate PAA 3%

LEVEL 6
480 min

Mixed Solution
Perform sterile concentrate Oxy 2%

LEVEL 6
480 min

Mixed Solution
Mucocit®-T branded mixture

LEVEL 6
480 min

50-01-1
Guanidine Hydrochloride

LEVEL 6
480 min

10127-02-3
Acridine orange

LEVEL 6
480 min

| | |
|--|--------------------|
| 77-86-1 Tris (hydroxymethyl) aminomethane Sat. solution | LEVEL 6 480 min |
| 127-09-3 Sodium acetate Sat. solution | LEVEL 6 480 min |
| 7664-93-9 Sulphuric Acid 50% | LEVEL 6 480 min |
| 76-03-9 Trichloroacetic acid 10% | LEVEL 6 480 min |
| 7803-57-8 Hydrazine monohydrate 80% | LEVEL 6 480 min |
| 7664-93-9 Sulphuric Acid 10% | LEVEL 6 480 min |
| 7681-52-9 Sodium Hypochlorite 13% | LEVEL 6 480 min |
| 1310-58-3 Potassium Hydroxide 40% | LEVEL 6 480 min |
| 7722-84-1 Hydrogen Peroxide 12% | LEVEL 6 480 min |
| 111-30-8 Gluteraldehyde 2.5% | LEVEL 6 480 min |
| 77-92-9 Citric acid 30% | LEVEL 6 480 min |
| 1310-73-2 Sodium Hydroxide, 50% | LEVEL 6 480 min |
| 1310-73-2 Sodium Hydroxide 40% | LEVEL 6 480 min |
| 7664-38-2 Phosphoric acid, 85% | LEVEL 6 480 min |

| | |
|--|------------------------------|
| 7664-38-2 Phosphoric Acid, 30% | LEVEL 6 480 min |
| 540-84-1 Iso-Octane 99% | LEVEL 6 480 min |
| 7722-84-1 Hydrogen Peroxide 30% | LEVEL 6 480 min DR 22% |
| 111-30-8 Glutaraldehyde 25% | LEVEL 6 480 min |
| 50-00-0 Formaldehyde 37% | LEVEL 6 480 min |
| 50-00-0 Formaldehyde 10% | LEVEL 6 480 min |
| 1239-45-8 Ethidium Bromide 5% | LEVEL 6 480 min |
| 7803-57-8 Hydrazine monohydrate 98% | LEVEL 4 180 min |
| 54-11-5 Nicotine 98% | LEVEL 4 151 min |
| 7647-01-0 Hydrochloric Acid 37% | LEVEL 3 116 min |
| Mixed Solution Klercide Premier - WFI 70/30 sterile IPA | LEVEL 3 106 min |
| 75-12-7 Formamide 99% | LEVEL 3 99 min |
| 110-54-3 n-Hexane 95% | LEVEL 3 97 min |
| Mixed Solution Klercide 70/30 sterile IPA | LEVEL 3 76 min |

| | |
|--|---------------------------|
| 78-83-1 Isobutanol 99% | LEVEL 3 76 min |
| 67-63-0 Isopropanol 70% | LEVEL 3 61 min |
| 108-39-4 m-Cresol 98.5% | LEVEL 2 59 min |
| 109-66-0 n-Pentane 98% | LEVEL 2 56 min |
| 598-75-4 Secondary isoamyl alcohol 98% | LEVEL 2 55 min |
| 108-87-2 Methylcyclohexane 99,9% | LEVEL 2 55 min |
| 7697-37-2 Nitric Acid, 50% | LEVEL 2 49 min |
| 110-82-7 Cyclohexane | LEVEL 6 480 min |
| 79-06-1 Acrylamide 40% | LEVEL 6 480 min |
| 67-63-0 Isopropanol 100% | LEVEL 2 43 min |
| 67-68-5 Dimethyl Sulfoxide 99% (DMSO) | LEVEL 2 42 min |
| Mixed Solution Klercide Premier - WFI 60/40 sterile Alcohol | LEVEL 2 41 min |
| 121-44-8 Triethylamine 99% | LEVEL 2 36 min |
| 64-17-5 Ethanol 70% | LEVEL 2 31 min |

| | |
|--|--------------------------|
| 108-95-2 Phenol 50% | LEVEL 1 24 min |
| 64-17-5 Ethanol 99.8% | LEVEL 1 16 min |
| 62-53-3 Aniline 99,9% | LEVEL 1 14 min |
| 95-63-6 1,2,4- Trimethylbenzene 98% | LEVEL 1 13 min |
| 97-88-1 Butyl methacrylate 99,9% | LEVEL 1 11 min |
| 1634-04-4 Methyl Tert Butyl Esther (MTBE) | LEVEL 1 11 min |
| 108-67-8 1,3,5-Trimethylbenzene 98% | LEVEL 1 10 min |
| 7664-39-3 Hydrofluoric Acid 40% | LEVEL 1 10 min |
| 127-19-5 Dimethyl Acetamide 99% | LEVEL 0 8 min |
| 74-97-5 Bromochloromethane | LEVEL 3 66 min |
| 64-19-7 Acetic Acid 100% | LEVEL 0 6 min |
| 7664-39-3 Hydrogen Fluoride 48% | LEVEL 0 6 min |
| 7664-93-9 Sulphuric Acid 95%-98% | LEVEL 0 6 min |
| 108-94-1 Cyclohexanone 99% | LEVEL 0 6 min |

| | |
|---|-------------------|
| 7664-39-3 Hydrofluoric Acid 48% | LEVEL 0 6 min |
| 67-56-1 Methanol 99,9% | LEVEL 0 5 min |
| 79-10-7 Acrylic acid 99% | LEVEL 0 4 min |
| 68-12-2 Dimethyl Formamide 99% | LEVEL 0 4 min |
| 71-36-3 Butanol 100% | LEVEL 2 39 min |
| 37143-54-7 1-Methoxy-2-propylamine 95% | LEVEL 0 3 min |
| 64-18-6 Formic acid 98,5% | LEVEL 0 3 min |
| 1336-21-6 Ammonium Hydroxide 25% | LEVEL 1 28 min |
| 1330-20-7 Xylene 98,5% | LEVEL 0 2 min |
| 96-47-9 2-Methyltetrahydrofuran 99,9% | LEVEL 0 1 min |
| 110-86-1 Pyridine | LEVEL 0 1 min |
| 100-51-6 Benzyl Alcohol | LEVEL 0 9 min |
| 75-56-9 Propylene oxide 99% | LEVEL 0 0 min |
| 109-99-9 Tetrahydrofuran 99.9% | LEVEL 0 0 min |

100-42-5
Styrene 99.9%

LEVEL 0
0 min

75-09-2
Dichloromethane 99%

LEVEL 0
0 min

107-21-1
Ethylene Glycol

LEVEL 0
0 min

67-66-3
Chloroform 99,8%

LEVEL 0
0 min

107-13-1
Acrylonitrile 99%

LEVEL 0
0 min

